

Validation of a Commercial Insulated Isothermal PCR-based POCKIT Test for Rapid and Easy Detection of White Spot Syndrome Virus Infection in *Litopenaeus vannamei*

Yun-Long Tsai¹, Han-Ching Wang², Chu-Fang Lo³, Kathy Tang-Nelson⁴, Donald Lightner⁴, Bor-Rung Ou⁵, Ai-Ling Hour⁶, Chuan-Fu Tsai¹, Cheng-Chi Yen¹, Hsiao-Fen Grace Chang¹, Ping-Hua Teng¹, Pei-Yu Lee¹*

1 Department of Research and Development, GeneReach Biotechnology Corporation, Taichung, Taiwan, 2 Institute of Biotechnology, National Cheng Kung University, Tainan, Taiwan, 3 Institute of Bioinformatics and Biosignal Transduction, National Cheng Kung University, Tainan, Taiwan, 4 Department of Veterinary Science and Microbiology, University of Arizona, Tucson, Arizona, United States of America, 5 Department of Animal Science and Biotechnology, Tunghai University, Taichung, Taiwan, 6 Department of Life Science, Fu-Jen Catholic University, Taipei, Taiwan

Abstract

Timely pond-side detection of white spot syndrome virus (WSSV) plays a critical role in the implementation of bio-security measures to help minimize economic losses caused by white spot syndrome disease, an important threat to shrimp aquaculture industry worldwide. A portable device, namely POCKITTM, became available recently to complete fluorescent probe-based insulated isothermal PCR (iiPCR), and automatic data detection and interpretation within one hour. Taking advantage of this platform, the IQ PlusTM WSSV Kit with POCKIT system was established to allow simple and easy WSSV detection for on-site users. The assay was first evaluated for its analytical sensitivity and specificity performance. The 95% limit of detection (LOD) of the assay was 17 copies of WSSV genomic DNA per reaction (95% confidence interval [CI], 13 to 24 copies per reaction). The established assay has detection sensitivity similar to that of OIE-registered IQ2000TM WSSV Detection and Protection System with serial dilutions of WSSV-positive *Litopenaeus vannamei* DNA. No cross-reaction signals were generated from infectious hypodermal and haematopoietic necrosis virus (IHHNV), *monodon* baculovirus (MBV), and hepatopancreatic parvovirus (HPV) positive samples. Accuracy analysis using700 *L. vannamei* of known WSSV infection status shows that the established assayhassensitivity93.5% (95% CI: 90.61–95.56%) and specificity 97% (95% CI: 94.31–98.50%). Furthermore, no discrepancy was found between the two assays when 100 random *L. vannamei* samples were tested in parallel. Finally, excellent correlation was observed among test results of three batches of reagents with 64 samples analyzed in three different laboratories. Working in a portable device, IQ PlusTM WSSV Kit with POCKIT system allows reliable, sensitive and specific on-site detection of WSSV in *L. vannamei*.

Citation: Tsai Y-L, Wang H-C, Lo C-F, Tang-Nelson K, Lightner D, et al. (2014) Validation of a Commercial Insulated Isothermal PCR-based POCKIT Test for Rapid and Easy Detection of White Spot Syndrome Virus Infection in *Litopenaeus vannamei*. PLoS ONE 9(3): e90545. doi:10.1371/journal.pone.0090545

Editor: Baochuan Lin, Naval Research Laboratory, United States of America

Received October 4, 2013; Accepted February 3, 2014; Published March 13, 2014

Copyright: © 2014 Tsai et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Funding: This work was supported in part by a grant from Central Taiwan Science Park, National Science Council, Taiwan, Grant no. 302403501 (http://www.ctspb. most.gov.tw/english/00home/home.aspx?v=20). The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript

Competing Interests: YLT, CFT, CCY, HFGC, PHT, and PYL are employed by GeneReach Biotechnology Corporation. However, this does not alter the authors' adherence to all the PLOS ONE policies on sharing data and materials.

* E-mail: pylee@genereachbiotech.com

Introduction

White spot syndrome, an OIE-listed disease, has resulted in severe economic losses in the shrimp aquaculture industry worldwide [1]. The disease is caused by white spot syndrome virus (WSSV), an enveloped large double-stranded DNA virus belonging to genus *Whispovirus* in the family *Nimaviridae*. Late post-larvae, juvenile and adult penaeid shrimp could all be infected by WSSV both vertically and horizontally. WSSV replicates preferentially in the cuticular epithelium and subcuticular connective tissues, as well as in the antennal gland, and haematopoietic tissues in shrimp. Factors such as rapid salinity and temperature changes can help activate WSSV in apparently healthy carriers and trigger the onset of the disease. Once activated, the disease could progress quickly and causes 100% mortality within 5 days [2,3]. In addition, WSSV could be found in a broad range of freshwater and marine crustaceans, including shrimp, crabs, crayfish and

lobsters [4,5,6]. Pre-screening of WSSV-free broodstock or larvae and regular surveillance of WSSV infection are important strategies to prevent WSSV infection outbreaks. Samples from pleopods, gills, haemolymph, stomach or abdominal muscle are recommended for diagnostic testing in shrimp [2].

Critical biosafety measurements could be implemented effectively when pond-side WSSV detection could be carried out to provide rapid, sensitive and specific test results. Polymerase chain reaction (PCR)-based assays, including basic PCR, nested PCR, and real-time PCR [7] assays, have been developed and demonstrated to be highly sensitive and specific for WSSV detection. However, these assays generally are not suitable for onsite applications because they require a trained technician and a relatively complicated thermocycler to operate. Furthermore, risks of cross contamination are relatively high due to the requirement of post-amplification processing of the reaction to detect PCR products for all PCR assays except for real-time PCR, which

Table 1. Assessment of analytical sensitivity of IQ Plus[™] WSSV Kit with POCKIT system.

Copies/reaction	Template					
	pWSSV1*	WSSV genomic DNA				
10 ³	100.00% (16/16)	100% (16/16)				
10 ²	100.00% (16/16)	100% (16/16)				
10 ¹	43.75% (7/16)	62.50% (10/16)				
10°	18.75% (3/16)	25.00% (4/16)				
10 ⁻¹	0.00% (0/16)	0.00% (0/16)				

*pWSSV1, standard WSSV plasmid. doi:10.1371/journal.pone.0090545.t001

nevertheless requires a relatively expensive device for signal detection and processing.

Recently, insulated isothermal PCR (iiPCR) was developed successfully to work in a simple thermally baffled device [8]. Within a specially designed cylindrical vessel (R-tubeTM; GeneReach Biotech, Taichung, Taiwan), iiPCR works by cycling reactions through temperature gradients established by the Rayleigh-Bénard convection when the tube is heated at the bottom with a single heating source. The three PCR steps, namely denaturation, annealing and extension, can be completed at different zones within the vessel. A WSSV iiPCR assay was reported to successfully produce detectable amplicons in 30 minutes with notable sensitivity and specificity [9]. However, the requirement of gel analysis to detect amplicons led to high risks of the assay amplicon cross-contamination, which was eliminated subsequently by incorporating hydrolysis probe technology into iiPCR and integration of an optical detection module into the iiPCR device [10]. An iiPCR device, namely POCKITTM, is now commercially available to perform the fluorescent probe-based iiPCR and provide test results automatically.

Taking advantage of the POCKIT system, IQ PlusTM WSSV Kit with POCKIT system was developed for on-site detection of WSSV in shrimp samples. Following procedures described in OIE Validation and Certification of Diagnosis Assays [11], the assay was evaluated and validated for WSSV detection in *Litopenaeus vannamei*. We report results for the assessment of analytical sensitivity, analytical specificity, repeatability, diagnostic accuracy, and proficiency of the established assay.

Materials and Methods

Sample Collection

L. vannamei samples were purchased from shrimp farms in Taiwan. All shrimp were shipped in oxygen-filled bags to GeneReach's lab and maintained for at least one day before screening. For experimental infection, a WSSV isolate obtained from naturally infected Penaeus monodon in Taiwan in 1994 was used. The status of WSSV infection in L. vannamei was defined by both an OIE reference method [5] and OIE-registered WSSV diagnostic kit, IQ2000TM WSSV Detection and Protection System (DPS). Infected shrimp were stored at -20°C before further analysis.

For proficiency evaluation of IQ PlusTM WSSV Kit with POCKIT system, abdominal muscles of *L. vannamei* experimentally infected with WSSV were aliquoted and preserved in 95% ethanol. A total of 31WSSV-negative and 34 WSSV-positive specimens were selected by IQ2000TM WSSV DPS. Aliquots of these samples were mailed to the laboratories of Dr. Chu-Fang Lo (National Taiwan University, Taiwan), Dr. Bor-Rung Ou (Tung-

Hai University, Taiwan), and Dr. Donald V. Lightner (University of Arizona, Tucson, USA) and tested by two operators at each site.

Nucleic Acid Extraction

DNA was extracted by using IQ PlusTM Extraction Kit (GeneReach Biotech) as described in the user manual. Briefly, tissues (50 mg) of ectodermal and mesodermal origin (e.g. abdominal muscle, pleopods, or periopods) were ground thoroughly in 500 µl of Solution 1. After the addition of Solution 2 (500 µl), the mixture was centrifuged at maximum speed using a cubee TM (GeneReach Biotech) at room temperature for 1 min. Supernatant (500 µl) was added to the spin column and centrifuged at maximum speed at room temperature for 1 min. The spin column was washed once with 500 µl of Solution 2 before DNA was eluted with 200 µl of Solution 3.

Preparation of WSSV Genomic DNA

Purified WSSV virions were prepared as described [12]. The concentration of purified DNA was determined by UV spectrophotometry and calculated on the basis of published genome size of 307,287 bp (GenBank accession no. AF440570).

IQ PlusTM WSSV Kit with POCKITsystem

The reaction in POCKITTM, an iiPCR-compatible instrument, completed a run in about one hour with the test readouts displayed on the monitor. Probe signals were collected during the reactions without the need to open the reaction vessels. POCKITTM was developed to include two separate channels, i.e. wavelengths 520 nm and 550 nm, on the basis of the optical detection module described previously for the original iiPCR device [8,10]. Signals passing through the 520- and 550-nm channels were collected by an integrated circuits controller-regulated CMOS sensor. Signalto-noise (S/N) ratios were calculated by dividing light signals collected after iiPCR by those from before iiPCR [10]. Limit of blank and positive cut-off for $POCKIT^{TM}$ were assigned on the basis of S/N ratios of numerous NTC (non-template control) and positive iiPCR reactions for various targets according to the reference method [13] and validated by analytical specificity and sensitivity tests (data now shown). Based on default thresholds of S/N ratios, results were converted automatically to "+", "-", or "?" and shown on the display screen. A "?" result indicated that the signals were ambiguous and the sample should be tested again.

Based on iiPCR and fluorescent probe hydrolysis for signal detection [10], IQ PlusTM WSSV Kit with POCKIT system (GeneReach Biotech) was designed to work in POCKITTM. Target-specific PCR primers and probe for WSSV were designed to interact with a segment of the *major capsid protein* gene of WSSV. To monitor DNA extraction and iiPCR reactions, primer pair and

T			WSSV-positive sample						~~~	~~~
1a	Target		10-1	10-2	10-3	10-4	10-5	NTC	SPF	STD
	signal nm)	3/3	3/3	3/3	3/3	3/3	0/3	0/3	0/3	3/3
	ignal nm)	3/3	3/3	3/3	3/3	3/3	3/3	0/3	3/3	3/3
M(bp)	10°	10-1	10-2	10-3 10-	4 10-5	101	STD 10 ²	103	NTC	SPF
					-	-,				
MARINE MARINE						-				
	-	M0000		100		MINIS				

Figure 1. Sensitivity analysis of IQ PlusTM WSSV Kit with POCKIT System. WSSV-positive *L. vannamei* DNAs were diluted in 200 ng/ml genomic DNA prepared from SPF *L. vannamei* and subjected to analysis with IQ PlusTM WSSV Kit with POCKIT System (A) and WSSV IQ2000 PDS (B). The signals from WSSV and IC amplification were collected by the 520-nm and 550-nm channels, respectively. (B) Positive standard DNA provided by the kit was serially diluted (10³, 10² and 10¹ copies) as suggested by the manufacturer. Nested PCR Amplicons generated by WSSV IQ2000 PDS were analyzed by agarose gel electrophoresis and ethidium bromide staining. The kit was designed to generate three amplicons (arrows) from WSSV genome and one from host shrimp DNA (arrow head). Number of target amplicon bands correlates positively with the initial concentrations of target DNA template. NTC, negative control (ddH2O); SPF, DNA of SPF *L. vannamei*; STD, positive control plasmid; M, DNA weight molecular markers. doi:10.1371/journal.pone.0090545.g001

probe targeting shrimp nuclear 18S rRNA gene were included to serve as the internal control (IC). All primers, probes, and dNTPs

were provided in a lyophilized Premix pellet. To assemble the reaction, the Premix pellet was first rehydrated in 50 µl of Premix




Figure 2. Specificity analysis of IQ PlusTM WSSV Kit with POCKIT System using WSSV, IHHNV, MBV, and HPV-positive samples. Genomic DNA extracts of diseased and healthy shrimp samples were subjected to analysis by IQ2000TM WSSV (A), IHHNV (B), MBV (C), and HPV (D) DPS assays, as well said PlusTM WSSV Kit with POCKIT system. Positive standard DNA provided by the kit was serially diluted (10³, 10² and 10¹ copies) as suggested by the manufacturer. Nested PCR amplicons were analyzed by agarose gel electrophoresis and ethidium bromide staining. The kits were designed to generate three amplicons (arrows) from target viral genome. The number of the amplicon bands positively correlates with the concentrations of the starting DNA template. The same set of samples was analyzed by IQ PlusTM WSSV Kit with POCKIT System analysis (E). M, DNA weight molecular markers (bp); N, negative control (ddH₂O); Std, standard plasmid. doi:10.1371/journal.pone.0090545.q002

4/4 Operator no. 2 Sample no. 0/4 0/4 0/4 0/4 0/4 0/4 2. Intra- and inter-assay repeatability of IQ PlusTM WSSV Kit with POCKIT system 4/4 4/4 4/4 4/4 4/4 Operator no. 1 Sample no. 0/4 0/4 0/4 0/4 0/4 0/4 ě ۲ Table Day

4/4 4/4 4/4 4/4 4/4

4/4

4/4

4 4

4/4 4/4

Buffer B. After 5 µl of extracted DNA were added to the Premix mixture, 50 µl of the final reaction mixture were transferred to an R-tubeTM (GeneReach Biotech). The R-tubeTM was sealed with a cap, spun for 10 sec in a cubee $^{\rm TM}$ (GeneReach Biotech), and placed into a POCKIT $^{\rm TM}$. The reaction was completed in one hour. Target and IC signals were collected through the 520- and 550-nm channels, respectively.

Nested PCR

IQ2000TM WSSV, IHHNV, MBV and HPV DPSs (GeneReach Biotech) were used to screen shrimp samples for the presence of target pathogens. Reactions were assembled and carried out as described in the manual. Briefly, the program for the first WSSV PCR included denaturation at 94°C for 2 min, 15 cycles of 94°C for 20 sec, 62°C for 20 sec and 72°C for 30 sec, and a cycle of 72°C for 30 sec and 20°C for 20 sec. The second PCR included 30 cycles of 94°C for 20 sec, 62°C for 20 sec and 72°C for 30 sec and one cycle of 72°C for 30 sec and 20°C for 20 sec. Nested PCR products were analyzed by electrophoresis on 2% agarose gels and stained with ethidium bromide.

Statistics

Statistical probit analysis, a non-linear regression model, was performed using commercial software SPSS 14.0 (SPSS Inc., Chicago, Illinois, USA) to determine limit of detection (LOD) with 95% confidence. Kappa statistic test (κ) was used to test the level of agreement between two different assays. Test results from different laboratories and operators in proficiency studies were analyzed by Chi-squire test and Fisher's exact test (for numbers less than 5) for homogeneity.

Results

IO PlusTM WSSV Kit with POCKIT System has Sensitivity Similar to Nested PCR

The recently developed POCKITTM device is designed specifically to provide optimal conditions for fluorescent probebased iiPCR. The device has the capacity to detect 520-nm and 550-nm signals simultaneously to detect multiplex signals. Taking advantage of the POCKIT system, a multiplex WSSV assay (IQ Plus TM WSSV Kit with POCKIT system) that could detect WSSV DNA and a shrimp IC target simultaneously was developed. Inclusion of an IC in the assay helps control for DNA extraction and subsequent enzymatic reactions.

First, to evaluate analytical sensitivity of IQ $Plus^{TM}$ WSSV Kit with POCKIT system, 10-fold serial dilutions of either linearized WSSV standard plasmid (pWSSV1) or WSSV genomic DNA of known copy numbers were analyzed. Analysis of 16 replicates of each dilution of standard plasmid shows that 100%, 43.75% and 18.75% of the 10^2 -, 10^1 -, and 10^0 -copyreactions, respectively, produced positive WSSV (520 nm) signals (Table 1). The LOD [95% hit rate] determined by probit regression analysis for standard plasmid was 24 copies per reaction (95% confidence interval [CI], 16 to 70 copies per reaction). Similarly, positive signals were generated from 100%, 62.5% and 25% of the reaction containing 10², 10¹ and 10⁰ copies of WSSV genomic DNA, respectively. The LOD determined by probit analysis with WSSV genomic DNA was17 copies per reaction (95% confidence interval [CI], 13 to 24 copies per reaction).

Furthermore, to compare the sensitivity of IQ PlusTM WSSV Kit with POCKIT system with that of IQ2000TM WSSV Detection and Prevention System (DPS), DNA extract prepared from a WSSV-positive L. vannamei was serially diluted using 10-fold increments in 20 ng/µl genomic DNA of SPF L. vannamei and

doi:10.1371/journal.pone.0090545.t002

Table 3. Accuracy test of IQ PlusTM WSSV Kit with POCKIT system using defined WSSV-positive and WSSV-negative *L. vannamei* samples.

	IQ2000 [™] WSSV DPS			
	Positive	Negative	Total	
Positive	374	9	383	
Negative	26	291	317	
Total	400	300	700	
	Negative	Positive Positive 374 Negative 26	Positive Negative Positive 374 9 Negative 26 291	

doi:10.1371/journal.pone.0090545.t003

subjected to analysis by both assays simultaneously. In iiPCR POCKIT assay, detection of the IC signals implied that both nucleic acid extraction and iiPCR steps were successful. Figure 1A shows that 550-nm IC signals were detected in all reactions containing shrimp DNA (Figure 1A, WSSV-positive sample and SPF). Among these, target WSSV signals (520 nm) were detected from 10^0 to 10^4 fold-diluted WSSV-positive samples. Detection limits of WSSV IQ2000 DPS were also reached at 10^{-4} dilution (Figure 1B), indicating that both assays had comparable detection limits to detect WSSV DNA in shrimp samples.

IQ PlusTM WSSV Kit with POCKIT System does not Cross-react with Important Shrimp DNA Viruses

In order to assess the specificity of the IQ $Plus^{TM}$ WSSV Kit with POCKIT system assay, L. vannamei infected strongly with important shrimp DNA viruses, including WSSV (sample No.1), IHHNV (infectious hypodermal and haematopoietic necrosis virus, sample No. 2), MBV (monodon baculovirus, sample No. 3), and HPV (hepatopancreatic parvovirus, sample No. 4), were selected by nested PCR-based IQ2000 assays (Figure 2A-2D). A healthy shrimp was also included (sample No. 5). The same set of samples was tested using IQ Plus WSSV Kit with POCKIT system (Figure 2E). IC signals, but no WSSV signals, were detected from the healthy shrimp. WSSV signals were detected only from the WSSV-positive but not from IHHNV-, MBV- or HPVpositive samples, which were also all tested IC positive (550 nm), suggesting that IQ PlusTM WSSV Kit with POCKIT system did not cross-react with these shrimp DNA viruses nor with host genome.

Precision Test

To evaluate intra-assay and inter-assay repeatability of IQ PlusTM WSSV Kit with POCKIT system, DNA extracted from one WSSV-negative and two WSSV-positive *L. vannamei* samples were analyzed in four replicates by two operators on two days. Table 2 shows that 100% agreement was found within and between tests of both WSSV (520 nm) and IC (550 nm) for each sample, demonstrating high precision of IQ PlusTM WSSV Kit

with POCKIT system for detection of WSSV DNA in shrimp samples.

Accuracy Test

IQ PlusTM WSSV Kit with POCKIT system was compared with IQ2000TM WSSV DPS for its accuracy in detecting WSSV DNA in shrimp samples. Sensitivity and specificity were determined according to the agreement between the two assays using positive and negative samples, respectively. First, 400 positive (including at least 100 lightly infected specimens) and 300 negative reference animals identified previously by IQ2000TM WSSV DPS (data not shown) were tested in parallel by both IQ PlusTM WSSV Kit with POCKIT system and IQ2000TM WSSV DPS. Results are summarized in Table 3. Kappa statistics performed to generate 95% confidence intervals indicates sensitivity: 93.5% [95% confidence interval (CI): 90.61–95.56%] and specificity: 97.0% (95% CI: 94.31–98.50%).

Furthermore, 100 un-defined shrimps obtained randomly from a local farm were tested. Table 4 shows that 100% agreement was observed between IQ2000TM WSSV DPS and IQ PlusTM WSSV Kit with POCKIT system for both sensitivity and specificity, with 19% of the samples tested positive by both assays.

Proficiency Test

Proficiency tests that involved three laboratories were performed to assess inter-lot, inter-laboratory and inter-operator reproducibility. Aliquots of 64 shrimp samples and three different batches of IQ Plus TM WSSV Kit with POCKIT system were shipped to the laboratories and repeated independently by two operators at each site. Six test results were obtained for each sample and analyzed by Chi-square test and Fisher's exact test for homogeneity evaluation of the kit. A total of 384 tests were done in each laboratory. Little heterogeneity was found among these three laboratories ($\chi^2 = 0.54$ and $\Lambda^2 = 0.76$, Chi-square test; $\Lambda^2 = 0.81$, Fisher's exact test). In addition, no significant difference was found among lots (Table 5) nor operators (Table 6), indicating the established kit performed with great reproducibility.

Table 4. Accuracy test of IQ PlusTM WSSV Kit with POCKIT system using undefined *L. vannamei* samples.

		IQ2000 [™] WSSV DPS			
		Positive	Negative	Total	
IQ Plus TM	Positive	19	0	19	
WSSV Kit with	Negative	0	81	81	
POCKIT system	Total	19	81	100	

doi:10.1371/journal.pone.0090545.t004

Table 5. Analysis of inter-lot reproducibility of IQ PlusTM WSSV Kit with POCKIT system.

Labs	Lots	DF	Total number	χ^2	P-value (χ^2)	P-value*
A	3	2	384	1.70	0.43	0.53
В	3	2	384	0.42	0.81	0.95
С	3	2	384	0.11	0.95	1

DF, degree of freedom. *Fisher's exact test. doi:10.1371/journal.pone.0090545.t005

Discussion

In this study, extensive evaluation and validation studies were performed to demonstrate the sensitivity, specificity, precision and accuracy of IQ PlusTM WSSV Kit with POCKIT system. Test results show that this system performed with great analytical sensitivity and specificity. The assay could detect as low as 17 copies of WSSV genomic DNA per reaction (95% LOD) and reach detection limits similar to those of the nested PCR-based IQ2000TM WSSV DPS, an OIE-registered WSSV detection assay. Furthermore, IQ PlusTM WSSV Kit with POCKIT system did not cross-react with other important shrimp DNA pathogens, including IHHNV, MBV, and HPV.

The POCKITTM device is developed to detect and measure fluorescent signals generated during PCR amplification to generate signal/noise (S/N) ratios. The S/N ratios are grouped into positive or negative readouts according to the built-in thresholds. Compared to OIE-registered $IQ2000^{TM}$ WSSV DPS, IQ PlusTM WSSV Kit with POCKIT system showed similar diagnostic accuracy on the basis of 700 prescreened L. vannamei with known WSSV infection status, i.e. sensitivity: 93.5% [95% CI: 90.61–95.56%] and specificity: 97.0% [95% CI: 94.31– 98.50%]) in this study. In addition, analysis of 100 randomly selected samples shows the established assay had 100% sensitivity and specificity. Therefore, positive and negative readouts derived from the built-in S/N threshold correlated well with the results from IQ2000TM WSSV DPS, supporting that the default S/N thresholds in POCKITTM for fluorescent probe-based iiPCR have been set properly.

This is the first report with rigorous validation studies to demonstrate that an assay based on the newly developed iiPCR POCKIT system could detect WSSV in shrimp sample with accuracy similar to conventional nested PCR. Compared to IQ2000TM WSSV DPS, different degrees of agreement were found in tests using 700 previously defined and 100 random samples for IQ PlusTM WSSV Kit with POCKIT System (Tables 3 and 4). Sensitivity of 93.5% and specificity of 97.0% were found from the analysis of 700 previously defined samples (Table 3), whereas100% agreement for both sensitivity and specificity was

observed with 100 randomshrimp samples (Table 4). The minor discrepancy found between the results of these two tests could be attributed to different sampling sizes and possible variations within sampling population, such as levels of infection. According to test results of IQ2000TM WSSV DPS, the previously defined and random populations were composed of different percentage of WSSV-negative samples (42.9% and 81%, respectively). Although validation tests were performed only with *L. vannamei* samples in this study, positive IC signals were also obtained from *P. monodon* (data not shown). Even though the 18S rRNA gene is conserved among crustaceans, whether the IC reaction in IQ PlusTM WSSV Kit with POCKIT system would work with a particular host should be evaluated before the assay is to be used for WSSV screening.

The IQ PlusTM WSSV Kit with POCKIT system requires simple steps and could be completed in one hour in a portable POCKITTM device without any post-amplification sample processing. Assembly of the reaction involves 3 steps, namely rehydration of the lyophilized reagents, addition of sample nucleic acids, and transfer of reaction mixture into reaction vessels. Automatic readouts generated with the default algorithm of POCKITTM are shown on the screen immediately after reaction completes, eliminating any manual data processing step to allow the assay to be performed by anyone with basic protocol training. In combination with a quick nucleic acid extraction kit, it takes less than 90 min from sample to result. However, the IQ PlusTM WSSV Kit with POCKIT system does require specialized reagents and equipment, including a mini-centrifuge and the POCKITTM device, to perform.

A variety of PCR-based assays have been developed for WSSV detection [7,14,15,16,17,18,19,20,21,22]. The requirements of either post-amplification process (conventional PCR) or sophisticated and expensive machine (real-time PCR) have limited their application to only laboratories with adequate equipment and trained technicians. In order to meet the requirements of a point-of-need WSSV detection assay, an iiPCR assay based on fluorescent probe hydrolysis was demonstrated to detect WSSV sensitively and specifically [10]. However, it was a singplex reaction performed in a prototype iiPCR device, which could

Table 6. Analysis of inter-operator reproducibility of IQ PlusTM WSSV Kit with POCKIT system in three laboratories.

Labs	Operators	DF	Total tests	χ^2	p-value (χ²)	p-value*
A	1,2	1	384	0	1	1
В	1,2	1	384	0.62	0.43	0.60
С	1,2	1	384	0.06	0.81	1

DF, degree of freedom. *Fisher's exact test.

doi:10.1371/journal.pone.0090545.t006

process signals through only one optical channel (520 nm) [10]. The recently launched POCKITTM included an optical system to detect 520-nm and 550-nmsignals separately to allow multiplex iiPCR, based on the prototype iiPCR device. The multiplex WSSV iiPCR described herein is the first assay to take full advantage of the POCKITTM device.

Suboptimal nucleic acid extraction and iiPCR conditions could lead to false-negative results. The IC reaction included in IQ PlusTM WSSV Kit with POCKIT system targets the shrimp's nuclear 18S rRNA gene, enabling users to monitor performance of the nucleic extraction and iiPCR reaction steps. Accordingly, the presence of IC signal and the absence of target signal indicate that WSSV is not detectable in the samples. The absence of IC signal alerts users that problems have occurred at DNA extraction and/or iiPCR stages, mitigating risks of obtaining false negative results.

Recently, various loop-mediated isothermal amplification (LAMP) assays have been developed for WSSV detection [23,24,25,26,27]. Requiring a simple incubator, LAMP also has the potential to be applied to point-of-need detection of aquaculture pathogens [28,29]. WSSV LAMP assays generally complete within one hour at around 63°C and have sensitivity limit of approximately 10 to 100 copies of template DNA. However, it is in general more challenging to optimize and multiplex LAMP assays because they require four to six primers to hybridize with six to eight areas in one target area. Detection of target amplicon-specific signals would reduce greatly the risks of false positive results derived from non-specific LAMP amplicons [23,26,30,31]. However, to our best knowledge, portable devices for target amplicon-specific detection are not available currently for field application of LAMP assays.

WSSV-like sequences, occupying around 20% of the *P. monodon* genome, has been found to be present throughout the shrimp genome in a genome sequencing study [32]. This is a cause for concern that the probability of false-positive diagnosis of white spot syndrome disease is likely to rise from cross reactivity of PCR primer and probe with target homologs within shrimp genome. Rigorous sequence BLAST analysis of the WSSV target sequence has been regularly carried out against sequences available in GenBank database, including all crustaceans. Furthermore, sequences of significant similarity with the target sequence have not been found in the *P. monodon* genomic sequences reported by Huang et al., 2011. Possible false positive detection of WSSV due to such sporadic DNA recombination events will be monitored continuously in the field, to alert us to implement proper adjustment of the assay in time.

High through-put real-time PCR assay or gel-based conventional nested PCR assays are recommended for laboratories staffed

References

- Sanchez-Paz A (2010) White spot syndrome virus: An overview on an emergent concern. Vet Res 41: 43.
- Lo CF, Ho CH, Chen CH, Liu KF, Chiu YL, et al. (1997) Detection and tissue tropism of white spot syndrome baculovirus WSBV) in captured brooders of *Penaeus monodon* with a special emphasis on reproductive organs. Dis Aquat Org 30: 53–72.
- Lotz JM, Soto MA (2002) Model of white spot syndrome virus (WSSV) epidemics in Litopenaeus vannamei. Dis Aquat Organ 50: 199–209.
- Hameed AS, Balasubramanian G, Musthaq SS, Yoganandhan K (2003) Experimental infection of twenty species of Indian marine crabs with white spot syndrome virus (WSSV). Dis Aquat Organ 57: 157–161.
- Lo CF, Ho CH, Peng SE, Chen CH, Hsu HC, et al. (1996) White spot syndrome baculovirus (WSBV) detected in cultured and captured shrimp, crabs and other arthroponds. Dis Aquat Org27: 215–225.
- Syed Musthaq S, Sudhakaran R, Balasubramanian G, Sahul Hameed AS (2006) Experimental transmission and tissue tropism of white spot syndrome virus (WSSV) in two species of lobsters, *Panulirus homarus* and *Panulirus omatus*. J Invertebr Pathol 93: 75–80.

with skilled technicians and required to process large numbers of samples for disease diagnosis and surveillance purposes. These methods could be afforded mostly by large-scale shrimp farms, and government-supported or commercial diagnostic laboratories. Shipping samples to be tested at centralized laboratories usually results in delayed responses and extra costs, especially for facilities at remote areas. The IQ Plus TM WSSV Kit with POCKIT system, a diagnosis assay allowing pond-side detection of WSSV, would help shrimp farmers and local offices to respond to disease outbreaks in an efficient and timely manner. For field users, equipment and accessories required to run the assay (POCKITTM, a mini-centrifuge, pipettes, and pipette tips) are combined into a mobile package (POCKITTM Xpress) to allow great mobility of the system. Compared to shipping samples to centralized laboratories for WSSV diagnosis, the POCKITTM assay could significantly lower the costs and shorten the sampling-to-result turn-around time from days to a few hours. Testing by IO PlusTM WSSV Kit with POCKIT system costs around US\$10 per sample, which is relatively inexpensive in comparison to the costs of sending samples to be diagnosed by standard and/or real-time PCR assays at a laboratory, which could cost more than US\$50 for each target pathogen plus fees for handling and shipping.

Allowing eight reactions per run, the field-deployable POCK-ITTM system is designed to fulfill the need of small-scale or branch facilities where, with a good sampling plan, high sample throughputs are not required. Furthermore, sensitivity of the POCKITTM system makes it possible to detect WSSV at early stages, enabling users to take appropriate bio-security measures in a timely manner. In addition, the POCKITTM system would facilitate on-site investigation of WSSV outbreaks, where relatively low numbers of known diseased animals are required to be tested.

Great reproducibility in WSSV detection by the IQ PlusTM WSSV Kit with POCKIT system was observed among the results generated by six operators at three laboratories, indicating excellent repeatability of the assay. This could be attributed partly to that the kit is available in a lyophilized format, which ensures stability of the reagents during the shipping and storage stages. Therefore, IQ PlusTM WSSV Kit assay has the potential to serve as a rapid, specific and sensitive tool for routine point-of-need detection of WSSV for shrimp farming industry.

Author Contributions

Conceived and designed the experiments: YLT HFGC PHT PYL. Performed the experiments: YLT CFL KT-N BRO CFT CCY. Analyzed the data: YLT HCW ALH PYL. Contributed reagents/materials/analysis tools: HCW. Wrote the paper: HCW CFL DL HFGC PYL.

- Durand SV, Lightner DV (2002) Quantitative realtime PCR for the measurement of white spot syndrome virus in shrimp. J Fish Dis 25: 381–389.
- Chang HFG, Tsai YL, Tsai CF, Lin CK, Lee PY, et al. (2012) A thermally baffled device for highly stabilized convective PCR. Biotech J 7: 662–666.
- Tsai YL, Lin YC, Chou PH, Teng PH, Lee PY (2012) Detection of white spot syndrome virus by polymerase chain reaction performed under insulated isothermal conditions. Journal Virological Methods 181: 134–137.
- Tsai YL, Wang HT, Chang HF, Tsai CF, Lin CK, et al. (2012) Development of TaqMan probe-based insulated isothermal PCR (iiPCR) for sensitive and specific on-site pathogen detection. PLoS One 7: e45278.
- OIE (2009) Principles of validation of diagnostic assays for infectious diseases.
 Manual of Diagnostic Tests for Aquatic Animals. Paris, France: Office international des épizooties.
- Xie X, Li H, Xu L, Yang F (2005) A simple and efficient method for purification of intact white spot syndrome virus (WSSV) viral particles. Virus Res 108: 63–
- Armbruster DA, Pry T (2008) Limit of blank, limit of detection and limit of quantitation. Clin Biochem Rev 29 Suppl 1: S49–52.

- 14. Kasornchandra J, Boonyaratpalin S, Itami T (1998) Detection of whitespot syndrome in cultured Penaeid shrimp in Asia: microscopic obserbation and polymerase chain reaction. Aquaculture 164: 243-251.
- 15. Lo CF, Leu JH, Ho CH, Chen CH, Peng SE, et al. (1996) Detection of baculovirus associated with white spot syndrome (WSBV) in penaeid shrimps using polymerase chain reaction. Dis of Aquat Org 25: 133-141.
- 16. Vaseeharan B, Jayakumar R, Ramasamy P (2003) PCR-based detection of white spot syndrome virus in cultured and captured crustaceans in India. Lett Appl Microbiol 37: 443-447.
- 17. Jian XF, Lu L, Chen YG, Chan SM, He JG (2005) Comparison of a novel in situ polymerase chain reaction (ISPCR) method to other methods for white spot syndrome virus (WSSV) detection in Penaeus vannamei. Dis Aquat Organ 67: 171-
- 18. Nunan LM, Lightner DV (2011) Optimized PCR assay for detection of white spot syndrome virus (WSSV). J Virol Methods 171: 318-321.
- 19. Belcher CR, Young PR (1998) Colourimetric PCR-based detection of monodon baculovirus in whole Penaeus monodon postlarvae. J Virol Methods 74: 21-29.
- 20. Kiatpathomchai W, Boonsaeng V, Tassanakajon A, Wongteerasupaya C, Jitrapakdee S, et al. (2001) A non-stop, single-tube, semi-nested PCR technique for grading the severity of white spot syndrome virus infections in Penaeus monodon. Dis Aquat Org 47: 235-239.
- 21. Natividad KD. Nomura N. Matsumura M (2008) Detection of White spot syndrome virus DNA in pond soil using a 2-step nested PCR. J Virol Methods
- Yuan L, Zhang X, Chang M, Jia C, Hemmingsen SM, et al. (2007) A new fluorescent quantitative PCR-based in vitro neutralization assay for white spot syndrome virus. J Virol Methods 146: 96–103. 23. Chou PH, Lin YC, Teng PH, Chen CL, Lee PY (2011) Real-time target-specific
- detection of loop-mediated isothermal amplification for white spot syndrome

- virus using fluorescence energy transfer-based probes. J Virol Methods 173: 67-
- 24. Kono T, Savan R, Sakai M, Itami T (2004) Detection of white spot syndrome virus in shrimp by loop-mediated isothermal amplification. J Virol Methods 115:
- 25. Mekata T, Sudhakaran R, Kono T, Supamattaya K, Linh NT, et al. (2009) Real-time quantitative loop-mediated isothermal amplification as a simple method for detecting white spot syndrome virus. Lett Appl Microbiol 48: 25-32.
- 26. Seetang-Nun Y, Jaroenram W, Sriurairatana S, Suebsing R, Kiatpathomchai W (2012) Visual detection of white spot syndrome virus using DNA-functionalized gold nanoparticles as probes combined with loop-mediated isothermal amplification. Mol Cell Probes 27: 71-79.
- Chaivisuthangkura P, Srisuk C, Rukpratanporn S, Longyant S, Sridulyakul P, et al. (2009) Rapid and sensitive detection of Penaeus monodon nucleopolyhedrovirus by loop-mediated isothermal amplification. J Virol Methods 162: 188-
- Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, et al. (2000) Loop-mediated isothermal amplification of DNA. Nucleic Acids Res 28: E63.
- Nagamine K, Hase T, Notomi T (2002) Accelerated reaction by loop-mediated isothermal amplification using loop primers. Mol Cell Probes 16: 223-229.
- 30. Mori Y, Hirano T, Notomi T (2006) Sequence specific visual detection of LAMP reactions by addition of cationic polymers. BMC Biotechnol 6: 3.
- Teng PH, Chen CL, Sung PF, Lee FC, Ou BR, et al. (2007) Specific detection of reverse transcription-loop-mediated isothermal amplification amplicons for Taura syndrome virus by colorimetric dot-blot hybridization. J Virol Methods 146: 317-326.
- 32. Huang SW, Lin YY, You EM, Liu TT, Shu HY, et al. (2011) Fosmid library end sequencing reveals a rarely known genome structure of marine shrimp Penaeus monodon. BMC Genomics 12: 242.